



HAPPY HANDS POWDER-FREE LATEX GLOVES – 4 MIL

Happy Hands Powder-Free Latex Gloves - 4 Mil

HH799 X-Small

HH800 Small

HH801 Medium

HH802 Large

HH803 X-Large

Federal Food, Drug, and Cosmetic Act (FD&C Act)



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Gloves](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tags: [gloves](#), [happy hands gloves](#), [latex gloves](#), [powder free gloves](#)

ADDITIONAL INFORMATION

Size 90, 89, 88, 91, 87

PRODUCT DESCRIPTION

Happy Hands Powder-Free Latex Gloves - 4 Mil

HH799 X-Small

HH800 Small



HH801 Medium

HH802 Large

HH803 X-Large

These premium gloves are made from natural rubber latex, offering formfitting comfort and fingertip dexterity. An offline manufacturing process gives the glove a low AQL (pinholes), and superiority over normal general purpose, B-grade gloves. In addition, Happy Hands latex gloves have a beaded cuff and reinforced finger thickness for excellent tear and puncture resistance.

Packaging: 100 gloves per dispenser box. 10 boxes per case.

Federal Food, Drug, and Cosmetic Act (FD&C Act)

The FD&C Act are laws passed by Congress giving authority to the U.S. Food and Drug Administration (FDA) to oversee the safety of food, drugs, and cosmetics. All Emerald Professional glove components are in compliance with provisions of the FD&C Act, which allow repeated use in direct contact with food.

ADDITIONAL INFORMATION

Size 90, 89, 88, 91, 87