



HOM1522 HEALTH-O-METER PEDIATRIC BALANCE BEAM SCALE WITH TRAY

HOM1522 Health-O-Meter Pediatric Balance Beam Scale with Tray

Weight Capacity: 130 lb / 65 kg

Heavy-duty Steel Base

Integrated Measuring Tape

2 Year Limited Warranty



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#), [Scales](#)

Tags: [beam scale](#), [pediatric scale](#)

1

PRODUCT DESCRIPTION

HOM1522 Health-O-Meter Pediatric Balance Beam Scale with Tray

Mechanical pediatric scales provide precision accuracy and ensure child safety. Durable steel bases and easy-to-clean trays guarantee long-lasting reliability. Child scales are ideal for pediatric and family practice physician offices.

Features:

- Heavy-duty Steel Base: Low-profile, non-tip steel base utilizes 3/32" steel with steel boxed levers and 3/16" hardened pivots and vee bearings for improved strength and durability
- Integrated Measuring Tape: Sanitary tray with built-in easy-to-read measuring tape allows user to measure both weight and height simultaneously



- Adjustable Unit Measurement Bar: Rotating kilogram/pound poise bars with engraved numerals allow easy reading that can withstand frequent use. Poise bar lock-in for error-free operation
- Factory calibrated to ensure accuracy
- Weight Capacity: 130 lb / 65 kg
- Graduation: 1 oz/20 g
- Measuring Tape: 0"-23"/0 cm-58 cm
- Tray Dimensions: 24 1/8" (w) x 14 1/2" (d) x 3 1/4" (h) / 613 mm (w) x 368 mm (d) x 83 mm (h)
- Product Dimensions: 20.75" x 24.12" x 22"
- Shipping Dimensions: 25.375" (l) x 20.625" (w) x 11.5" (h) / 3.48 Cube
- Shipping Weight: 31 lb
- 2 Year Limited Warranty