




HOM2842 HEALTH-O-METER LARGE PLATFORM DIGITAL BARIATRIC/ VETERINARY SCALE WITH REMOTE DISPLAY

HOM2842 Health-O-Meter Large Platform Digital Bariatric/ Veterinary Scale with Remote Display
Weight Capacity: 600 lb/270 kg
2 (Two) Year Limited Warranty



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#), [Scales](#)

Tags: [platform scale](#), [veterinary scale](#)

1

PRODUCT DESCRIPTION

HOM2842 Health-O-Meter Large Platform Digital Bariatric/ Veterinary Scale with Remote Display

DIGITAL VET SCALE

This scale will fit every need of the healthcare professional or veterinarian with its stainless steel platform and durable anti-skid mat. Ideal for healthcare facilities and veterinary offices, this scale can accommodate patients, animals, or medical linen carts up to 600 lb.

Features:

- Motion Detection: Scale has the ability to weigh accurately despite movement of the subject on the scale
- Durable Construction: Easy-to-clean stainless steel base with removable non-skid mat can



withstand frequent use

- Remote Display: Easy-to-read 7/8" LCD digital display provides versatility with the ability to be mounted on a wall or placed on a table (mounting bracket, 6 foot cable included)
- EMR Connectivity via Serial Port: Scale has the ability to interface with a computer, monitor, printer or other electronic device via Serial Port for fast, easy and accurate data transfer
- Weight Capacity: 600 lb/270 kg
- Graduation: .2 lb/.1 kg
- Large Low-Profile Platform: 22 1/4" (w) x 42" (d) x 2 1/2" (h)/565 mm (w) x 1067 mm (d) x 64 mm (h)
- Product Dimensions: 42" x 22.37" x 2.6"
- Unit Conversion: lb/kg
- Zero Out/Tare Functions
- Low Battery Indicator
- Dual Power Sources: Included 120V AC adapter or 6 C-cell batteries (not included)
- Shipping Dimensions: 47.5" (l) x 25.625" (w) x 5.25" (h)/3.70 Cube
- Shipping Weight: 61.75 lb
- Warranty: 2 (Two) Year Limited Warranty