



HOM553 DIGITAL PEDIATRIC/ VETERINARY SCALE

HOM553 Digital Pediatric/ Veterinary Scale
Weighs in kilograms and pounds/ounces
Weight capacity of 44lbs/20kg
Powered by 6 AA batteries



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#), [Scales](#)

Tags: [digital pediatric scale](#), [digital veterinary scale](#)

1

PRODUCT DESCRIPTION

HOM553 Digital Pediatric/ Veterinary Scale

High performance digital pediatric scale with motion-sensing technology to capture the weights of even the most active babies.

Weighs in kilograms and pounds/ounces

Weight capacity of 44lbs/20kg

Waterproof On/Off, Zero/Clear (tare), Hold/Release and Kg/lb buttons

Features include auto-off, auto-zero and low-battery indicator.

Tare memory

Measuring tape included (24in x 1/8in / 59cm x 1mm)

Powered by 6 AA batteries (included) or AC adapter (not included)

