




HYBRIDLX ROLLATOR TRANSPORT CHAIR

HybridLX Rollator Transport Chair

LX1000T Lumex Hybrid LX Rollator Transport Chair, Titanium

LX1000B Lumex Hybrid LX Rollator Transport Chair, Majestic Blue

Three (3) year limited warranty



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DURABLE MEDICAL](#), [Transport Chair](#)



Tags: [rollator](#), [transport chair](#)

ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Color	1542, 1541

VARIATIONS



Image	SKU	Price	Description	Color
				Titanium
				Majestic Blue

PRODUCT DESCRIPTION

HybridLX Rollator Transport Chair

LX1000T Lumex Hybrid LX Rollator Transport Chair, Titanium

LX1000B Lumex Hybrid LX Rollator Transport Chair, Majestic Blue

The HybridLX is an innovative rollator that combines the function of a wheeled walker and a transport chair into one lightweight design. Keeping the patient's desire for independence and comfort in mind, the HybridLX is designed with secure handgrips, a comfortable and wide seat and backrest, and the effortless design to transition from a rollator to a transport chair.

Other Features Include:

- Lightweight frame - durable, aluminum frame weighs only 22 lbs (this includes footrests) making it easier for users to move about freely
- Hybrid concept - designed to be used primarily as a rollator while easily transforming into a transport chair with removable footrests (adjustable in length) when user needs a short rest
- Ultra-Compact design - easy to assemble without the need for tools and folds effortlessly for convenient transport and storage
- Easy-to-use handbrakes - allow for short and longer period braking when used as rollator and transport chair (user must be accompanied when using as a transport chair)



- Flexible, soft polyester backrest and 18" wide mesh seat with pull up strap for easy folding provides comfort for user
- Large 8" wheels for easy rolling indoor and outdoor on rough surfaces
- Complete with storage pouch for personal items
- Available in attractive, contemporary colors: Titanium and Majestic Blue
- Designed to meet the requirements of HCPCS Codes: E0143 + E0156.
- Three (3) year limited warranty

ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Color	1542, 1541