



80071 CLINTON IMAGING TABLE FLAT TOP W/ DROP WINDOW

80071 Clinton Imaging Table Flat Top W/ Drop Window

Item No: 80071



SKU: N/A

Price:

Stock: Available on backorder

Categories: [IMAGING TABLES](#), [Ultrasound Table](#)

Tags: [imaging table](#), [ultrasound table](#)

1

PRODUCT DESCRIPTION

80071 Clinton Imaging Table Flat Top W/ Drop Window

Item No: 80071

CLINTON ONE PIECE FLAT TOP WITH DROP WINDOW; TRENDELENBURG/REVERSE TRENDELENBURG PLUS OR MINUS 15 ; EASY-TO-USE HAND CONTROL; ALL-WELDED BASE AND TOP FRAME; UNOBSTRUCTED BASE DESIGN; DUAL COLUMN LIFT SYSTEM ADDS RIGIDITY; PAPER DISPENSER INCLUDED; POWER SUPPLIES AND MOTORS MEET UL 60601-1 AND CAN/CSA C22.2 601.1; INDUSTRY LEADING 5-YEAR, LIMITED WARRANTY; 550 LBS LOAD CAPACITY UNDER NORMAL USE

1