



## LF1600 NON-POWERED STAND ASSIST

### LF1600 Non-Powered Stand Assist

Note: The Lumex Stand Assist is designed for residents who have the strength and stability to lift and support themselves. We recommend an electric sit-to-stand lift for residents with limited strength and stability.



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [Non-Powered Stand Assist](#), [PATIENT LIFTS](#)

**Tag:** [stand assist](#)

1

## PRODUCT DESCRIPTION

### LF1600 Non-Powered Stand Assist

#### LF1600 Lumex Stand Assist Patient Transport

The Lumex Stand Assist is a transport assistance unit that keeps the resident actively engaged in the transport process. Transport functions of all kinds are quick, convenient, and require minimal caregiver assistance. Users grasp the middle bar and use their own strength to pull themselves up into position. A padded split seat swings out to allow loading or unloading, then swings back to form a comfortable and stable transport seat. The Stand Assist is a convenient toileting alternative to a wheelchair.

- 400 lb safe working load
- Padded dual knee pads and seat pads for comfortable transport



- Stable base with locking rear casters
- Fits easily through narrow doorways
- Affordable alternative to electric sit-to-stand lifts

Note: The Lumex Stand Assist is designed for residents who have the strength and stability to lift and support themselves. We recommend an electric sit-to-stand lift for residents with limited strength and stability.

#### Specifications

Maximum Width.....25"

Maximum Height.....42.5"

Overall Length.....33.5"

Unit Weight.....61 lb

Front Base Height.....5"

**WARNING:** Personal injury - this device is designed for transport assistance of individuals between approximately 5'1" and 6'6" (155 cm and 198 cm) in height. Ensure the LF1600 properly fits the individual being transferred before use.