



LIBERTY HALF NO GAP BED RAIL

Liberty Half No Gap Bed Rail
 GF6590BH-1 BED RAILS LIBERTY HALF HIGH NO GAP 2 / PR
 GF6680B-1 Lumex Liberty Bed Rail Crossbars 2 / PR Gray



SKU: N/A
Price:
Stock: 10 in stock (can be backordered)
Categories: [Bed Rails](#), [Patient Room](#)
Tags: [bed rail](#), [no gap bed rail](#)

ADDITIONAL INFORMATION

Item No. 2 279, 278

VARIATIONS

Image	SKU	Price	Description	Item No. 2
				GF6680B-1



Image

SKU

Price

Description

Item No. 2



GF6590BH-1

PRODUCT DESCRIPTION

Liberty Half No Gap Bed Rail

GF6590BH-1 BED RAILS LIBERTY HALF HIGH NO GAP 2 / PR

GF6680B-1 Lumex Liberty Bed Rail Crossbars 2 / PR Gray

The Half No-Gap Head Rail, which provides a full head coverage of 31.5 inches, extends to eliminate the gap between the head board and the end of the rail. The non-telescoping tube overlap results in a very sturdy rail that can be used as a secondary handle and support. Designed to offer added protection, the Lumex Liberty Rails follow the FDA Hospital Bed System Dimensional Guidance to Reduce Entrapment and the Hospital Bed Safety Workgroup and FDA Dimensional Test Methods for Bed System Zones 1-4 to reduce bed entrapment.

Features:

Improved rail design decreases the spacing between the crossbars, thus reducing the risk of entrapment

New bed rail length closes open areas that can occur during full bed articulation

Rails now mount in fixed patent pending rail slots on Patriot™ bed deck to prevent sliding and moving — easy installation saves delivery and setup time and ensures proper rail placement

Spring-loaded, self-adjusting crossbraces feature safety tabs that prevent accidental disengagement when bumped or during application of side pressure

Patriot™ Liberty Bed Rails have been tested to the highest Medical Bed System standards

Designed to be used specifically with the Patriot™ and WhisperLite™ II full-electric beds that are commonly used as standard-height or low beds

Heavy-duty steel construction

Rails are powder-coated with an environmentally-friendly, long-lasting hammertone finish and

