



LIBERTY QUARTER LENGTH BED RAIL

Liberty Quarter Length Bed Rail

GF6580BH-1 Lumex Liberty Quarter Length Bed Rail 2 / PR Gray

GF6680B-1 Lumex Liberty Bed Rail Crossbars 2 / PR Gray



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Bed Rails](#), [Patient Room](#)

Tags: [bed rail](#), [liberty bed rail](#), [quarter bed rail](#)

ADDITIONAL INFORMATION

Item No. 2 280, 278

VARIATIONS

Image

SKU

Price

Description

Item No. 2



GF6680B-1



Image	SKU	Price	Description	Item No. 2
-------	-----	-------	-------------	------------



GF6580BH-1

PRODUCT DESCRIPTION

Liberty Quarter Length Bed Rail

GF6580BH-1 Lumex Liberty Quarter Length Bed Rail 2 / PR Gray

GF6680B-1 Lumex Liberty Bed Rail Crossbars 2 / PR Gray

The Quarter Rail, approximately 18 inches in length when mounted on the Patriot™ bed, can be used as a head assist, head rail, or foot rail. This rail provides minimal rail coverage and gives assistance when entering or exiting the bed. The Quarter Rail's non-telescoping inner and outer tubes overlap in conjunction with the Patriot™ mounting system, resulting in a very sturdy rail. Designed to offer added protection, the Lumex Liberty Rails follow the FDA Hospital Bed System Dimensional Guidance to Reduce Entrapment and the Hospital Bed Safety Workgroup and FDA Dimensional Test Methods for Bed System Zones 1-4 to reduce bed entrapment.

Features:

Improved rail design decreases the spacing between the crossbars, thus reducing the risk of entrapment

New bed rail length closes open areas that can occur during full bed articulation

Rails now mount in fixed patent pending rail slots on Patriot™ bed deck to prevent sliding and moving — easy installation saves delivery and setup time and ensures proper rail placement

Spring-loaded, self-adjusting crossbraces feature safety tabs that prevent accidental disengagement when bumped or during application of side pressure

Patriot™ Liberty Bed Rails have been tested to the highest Medical Bed System standards

Designed to be used specifically with the Patriot™ and WhisperLite™ II full-electric beds that are commonly used as standard-height or low beds

Heavy-duty steel construction



Rails are powder-coated with an environmentally-friendly, long-lasting hammertone finish and secure gripping surface

Tool-free installation

The Quarter Rail can be mounted on the bed's foot and the GF6590B-1 mounted on the bed's head to create the Patriot™ split rail system.

3 Year Limited Warranty

Designed to meet the requirements of HCPCS E0310

Specifications

Overall length: (full extension) 18"

Overall height: 13.5"

Product Weight: 12.75 lb Specifications are subject to change without notice.

The information provided herein regarding recommended usage is in fact a recommendation. Please consult your healthcare professional for proper guidance and selection of homecare bed rails to best suit intended user needs.

ADDITIONAL INFORMATION

Item No. 2 280, 278