




LUMEX® BARIATRIC EASY LIFT STS

Lumex® Bariatric Easy Lift STS
LF2090 LIFT SIT TO STAND 600 LB
DPL650BAT BATTERY LF1050 1090 2020 2090
DPL650PEN HAND CTL LF1050 1090 2020 2090
DPL650ADP CHARGER LF1050 1090 2020 2090



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Battery Powered Lifts](#), [PATIENT LIFTS](#)

Tags: [bariatric lift](#), [easy lift](#), [patient lift](#)



ADDITIONAL INFORMATION

Item No. 2 367, 365, 366, 364



VARIATIONS

Image	SKU	Price	Description	Item No. 2
				LF2090

PRODUCT DESCRIPTION

Lumex® Bariatric Easy Lift STS
 LF2090 LIFT SIT TO STAND 600 LB
 DPL650BAT BATTERY LF1050 1090 2020 2090
 DPL650PEN HAND CTL LF1050 1090 2020 2090
 DPL650ADP CHARGER LF1050 1090 2020 2090
 Heavy-gauge steel construction
 Designed to provide quick, safe assistance to those who have trouble standing
 Unique, ergonomic foot pedal base opening
 Large, easy to grip handles for easy maneuvering
 Lower base height maximizes stability and fits under most beds
 Rubber base coating protects furniture and walls
 Adjustable knee pad unit
 Quick release removable foot platform
 Optional buttock strap for patients requiring additional support
 24V DC motor allows quiet, smooth operation
 Ergonomically designed hand pendant with pendant hook
 Warning beep for low battery charge and weight capacity levels
 Easy access emergency stop button
 White powder-coated finish
 Accessory storage bag included for slings and charger
 Interchangeable battery works with other Lumex Easy Lifts
 Easy access emergency stop button



Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED

Designed to meet the requirements of HCPCS Code E0635

Warranty:

Frame and spreader bar three-year limited

Actuator two-year limited

Batteries one-year limited

ADDITIONAL INFORMATION

Item No. 2 367, 365, 366, 364