



LUMEX® EASY LIFT PATIENT LIFTING SYSTEM - BARIATRIC

Lumex® Easy Lift Patient Lifting System - Bariatric
LF1090 600 lb weight capacity
DPL650BAT BATTERY LF1050 1090 2020 2090
DPL650PEN HAND CTL LF1050 1090 2020 2090
DPL650ADP CHARGER LF1050 1090 2020 2090



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Battery Powered Lifts](#), [PATIENT LIFTS](#)

Tags: [battery-powered lift](#), [easy lift](#), [patient lift](#)

ADDITIONAL INFORMATION

Item No. 2 367, 365, 366, 372

VARIATIONS



Image	SKU	Price	Description	Item No. 2
				LF1090

PRODUCT DESCRIPTION

Lumex® Easy Lift Patient Lifting System - Bariatric

LF1090 600 lb weight capacity

DPL650BAT BATTERY LF1050 1090 2020 2090

DPL650PEN HAND CTL LF1050 1090 2020 2090

DPL650ADP CHARGER LF1050 1090 2020 2090

Heavy-gauge steel construction

Six-Point spreader bar with 360° rotation tilts to enhance comfort and safety

Uses Lumex six-point, four-point, or two-point slings

Unique, ergonomic foot pedal leg opening

Large, easy to grip handles for easy maneuvering

Lower base height maximizes stability and fits under most beds

Rubber base coating protects furniture and walls

24V DC motor allows quiet smooth operation

Ergonomically designed hand pendant with pendant hook

Warning beep for low battery charge and weight capacity levels

Easy access emergency stop button

Emergency manual lowering device

White powder-coated finish

Accessory storage bag included for slings and charger

Optional digital scale available

Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED

Designed to meet the requirements of HCPCS Code E0635

Warranty:

Frame and spreader bar three-year limited



Actuator two year limited
Batteries one year limited

ADDITIONAL INFORMATION

Item No. 2 367, 365, 366, 372