



LUMEX OVERBED TABLE

Lumex Overbed Table

GF9000-LSTMO Overbed Table - Single Top, Laminate, Medium Oak

GF9000-LSTMC Overbed Table - Single Top, Laminate

GF9000-LSTHM Overbed Table - Single Top, Laminate

GF9000-LSTCL Overbed Table - Single Top, Laminate

GF9000-TSTMO Overbed Table - Single Top, Thermofoil

GF9000-TSTMC Overbed Table - Single Top, Thermofoil

GF9000-TSTHM Overbed Table - Single Top, Thermofoil

F9000-TSTCL Overbed Table - Single Top, Thermofoil

GF9000-TDTMO Overbed Table - Flip Top, Thermofoil

GF9000-TDTMC Overbed Table - Flip Top, Thermofoil

GF9000-TDTHM Overbed Table - Flip Top, Thermofoil

GF9000-TDTCL Overbed Table - Flip Top, Thermofoil

GF9000-BSK Storage Basket - 10X6x9 With Adapter

GF9000-VNTY Vanity Drawer With Mirror

GF9000-UH Urinal Hook With Adapter

GF9000-CST Non-Locking Caster



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Overbed Table](#), [Patient Room](#)

Tags: [lumex table](#), [overbed table](#)

ADDITIONAL INFORMATION



Weight	N/A
Dimensions	N/A
Item No. 2	299, 304, 307, 295, 294, 293, 292, 303, 302, 301, 300, 298, 297, 296, 306, 305

VARIATIONS

Image	SKU	Price	Description	Item No. 2
				GF9000-LSTMO
				GF9000-LSTMC
				GF9000-LSTHM



Image

SKU Price Description

Item No. 2



GF9000-LSTCL



GF9000-TSTMO



GF9000-TSTMC



GF9000-TSTHM



Image

SKU Price Description

Item No. 2



F9000-TSTCL



GF9000-TDTMO



GF9000-TDTMC



GF9000-TDTHM



Image

SKU Price Description

Item No. 2



GF9000-TDTCL



GF9000-BSK



GF9000-VNTY



GF9000-UH



Image

SKU Price Description

Item No. 2



GF9000-CST

PRODUCT DESCRIPTION

Lumex Overbed Table

- GF9000-LSTMO Overbed Table - Single Top, Laminate, Medium Oak, Graphite Grey U Base
- GF9000-LSTMC Overbed Table - Single Top, Laminate, Medium Cherry Graphite Grey U Base
- GF9000-LSTHM Overbed Table - Single Top, Laminate, Hardrock Maple, Graphite Grey U Base
- GF9000-LSTCL Overbed Table - Single Top, Laminate, Candlelight, Graphite Grey U Base
- GF9000-TSTMO Overbed Table - Single Top, Thermofoil, Medium Oak, Graphite Grey U-Base
- GF9000-TSTMC Overbed Table - Single Top, Thermofoil, Medium Cherry Graphite Grey U-Base
- GF9000-TSTHM Overbed Table - Single Top, Thermofoil, Hardrock Maple, Graphite Grey U-Base
- GF9000-TSTCL Overbed Table - Single Top, Thermofoil, Candelight, Graphite Grey U-Base
- GF9000-TDTMO Overbed Table - Flip Top, Thermofoil, Medium Oak, Graphite Grey U-Base
- GF9000-TDTMC Overbed Table - Flip Top, Thermofoil, Medium Cherry Graphite Grey U-Base
- GF9000-TDTHM Overbed Table - Flip Top, Thermofoil, Hardrock Maple, Graphite Grey U-Base
- GF9000-TDTCL Overbed Table - Flip Top, Thermofoil, Candelight, Graphite Grey U-Base
- GF9000-BSK Storage Basket - 10X6x9 With Adapter
- GF9000-VNTY Vanity Drawer With Mirror
- GF9000-UH Urinal Hook With Adapter
- GF9000-CST Non-Locking Caster

The Lumex Overbed Table combines an attractive, ergonomic design with features that reduce infection control issues by minimizing areas in which dust and dirt can accumulate. Manufactured with the lowest base in the industry, our Overbed tables move easily under any hospital bed or recliner.



Features:

- Lowest profile base on the market designed to fit under any hospital bed
- Exposed steel base with no plastic cover allows for easy cleaning of base and casters
- Laminate or Thermofoil top that provides durability
- Available in 4 attractive finishes: Candlelight, Hardrock Maple, Medium Cherry, and Medium Oak
- Seamless Spill guard top holds up to 36 oz
- Easy, one-hand activated secondary surface reduces infection control issues and eliminates hidden spaces between surfaces for bacteria to collect
- Integrated cup holders can hold soda cans, water bottles and medication
- Optional Vanity drawer with mirror slides out from either side of the table
- Integrated accessory rail on both sides of column to easily transfer accessories such as cord wraps, IV hooks, urinal hooks and storage baskets from headwall rails
- Shock-absorbent gas spring cylinder for easy height adjustment
- Designed to support up to 100 lb of evenly distributed weight
- Two (2) Year Limited Warranty

Standard Dimensions

Single Table Top:

32" L x 17" W (81.3 cm x 43.2 cm)

Flip Table Top:

13.5" L x 17" W (34.3 cm x 43.2 cm)

Storage Table Top:

41.5" L x 19" W (105.4 cm x 48.3 cm)

Height Range:

28" Low - 44" High (71 cm - 111.8 cm)

ADDITIONAL INFORMATION

Weight N/A

Dimensions N/A

Item No. 2 299, 304, 307, 295, 294, 293, 292, 303, 302, 301, 300, 298, 297, 296, 306, 305