




MOIST HEAT THERAPY PACKS

Moist Heat Therapy Packs
2086 Standard Pack
2087 Half Size Pack
2088 Neck Pack



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Personal Care/Daily Living Aids](#), [Therapy & Rehab](#)




Tags: [heat therapy packs](#), [moist therapy packs](#), [therapy packs](#)

ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Item No. 4	484, 485, 486

VARIATIONS



Image	SKU	Price	Description	Item No. 4
				2086
				2087
				2088

PRODUCT DESCRIPTION

Moist Heat Therapy Packs

2086 Standard Pack 10" x 12"

2087 Half Size Pack 5" x 12"

2088 Neck Pack 9" x 24"

Moist heat pack is a simple and effective method of applying moist heat to a local area

Gives at least 30 minutes of therapeutic heat

Use with either a moist pack cover or a heavy towel



ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Item No. 4	484, 485, 486