



ORANGE MAX-GRIP TREAD RIBBED NITRILE GLOVES

Orange Max-Grip Tread Ribbed Nitrile Gloves


O3999 Medium

O4000 Large

O4001 X-Large

O4002 XX-Large

Federal Food, Drug, and Cosmetic Act (FD&C Act)



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Gloves](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tags: [gloves](#), [Max-Grip gloves](#), [nitrile gloves](#), [ribbed gloves](#)

ADDITIONAL INFORMATION

Size	90, 89, 91, 92
-------------	----------------

PRODUCT DESCRIPTION

Orange Max-Grip Tread Ribbed Nitrile Gloves

O3999 Medium

O4000 Large

O4001 X-Large

O4002 XX-Large



The Emerald Orange Max-Grip powder-free nitrile glove offers an innovative design giving tenacious grip with superior tactile feel. This latex-free 7 mil disposable glove has a beaded cuff and a 4 mil textured ridge to channel fluids away when gripping wet surfaces. The vibrant orange color makes this the ideal glove for PPE applications requiring high-vis clothing. Perfect for mechanics or anyone needing protection and security handling slick equipment. 50 gloves per dispenser box, 10 boxes per case.

Federal Food, Drug, and Cosmetic Act (FD&C Act)

The FD&C Act are laws passed by Congress giving authority to the U.S. Food and Drug Administration (FDA) to oversee the safety of food, drugs, and cosmetics. All Emerald Professional glove components are in compliance with provisions of the FD&C Act, which allow repeated use in direct contact with food.

ADDITIONAL INFORMATION

Size	90, 89, 91, 92
-------------	----------------