



PANASCOPE™ STETHOSCOPES-LIGHTWEIGHT

Panascope™ Stethoscopes-Lightweight
500 With Adult Chestpiece Black
500GY With Adult Chestpiece Grey
510 With Pediatric Chestpiece Black
510GY With Pediatric Chestpiece Grey



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DIAGNOSTIC INSTRUMENTS](#), [Stethoscopes](#)

Tags: [panascope stethoscopes](#), [stethoscopes](#)

ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Color	107, 110
size age	144, 145

VARIATIONS



Image	SKU	Price	Description	Color	size age
				Black	adult
				Black	child
				Grey	adult
				Grey	child

PRODUCT DESCRIPTION

Panascope™ Stethoscopes-Lightweight
500 With Adult Chestpiece Black



500GY With Adult Chestpiece Grey

510 With Pediatric Chestpiece Black

510GY With Pediatric Chestpiece Grey

Lightweight deluxe stethoscope offers an aluminum dual head chestpiece with a chrome-plated rotating stem.

Non-chill ring

22" latex-free tubing.

Adjustable binaurals are chrome-plated brass with an internal spring. Comes complete with additional diaphragm and eartips.

Lifetime Warranty.

ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Color	107, 110
size age	144, 145