



PATRIOT HOMECARE BEDS FULL-ELECTRIC/LOW BEDS

Patriot Homecare Beds Full-Electric/Low Beds

Total assembled bed weighs only 178 pounds (bed and motor).

Powerful DC motor for effortless height adjustment.

Hi-Lo motor comes pre-installed on foot section for easy set-up and disassembly.

Designed to meet the requirements of applicable HCPCS codes.

Maximum patient weight 350 lbs.

450 lb evenly distributed total weight capacity including patient, mattress, rails, and accessories.

Height Adjustment:

Casters in Low Position:

Min. 9.5", Max. 20" (at deck level)

Casters in Standard Position:

Min. 13", Max. 23.5" (at deck level)

Overall Bed Size: 87" L x 36" W.

Lift Range: Head - 60°, Knee - 30°.

No external pinch points!



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [BED](#), [Full Electric Beds](#)

Tags: [full-electric bed](#), [Homecare bed](#), [patient bed](#)

ADDITIONAL INFORMATION



Item No. 2

1719, 1737, 1720, 1721, 1722, 1723,
1724, 1725, 1726, 1727, 1728, 1729,
1730, 1731, 1732, 1733, 1734, 1735,
1736

VARIATIONS

Image	SKU	Price	Description	Item No. 2
				US0458
				US0458-RFPKG
				US0458-RPKG



Image	SKU	Price	Description	Item No. 2
				US0458PL
				US0458PL-220V
				US0458PL-RFPKG
				US0458PL-RFPKGHA



Image	SKU	Price	Description	Item No. 2
				US0458PL-RPKG
				US0458PL-RPKG HRA
				US0468-PKG
				US0468-RFPKG



Image	SKU	Price	Description	Item No. 2
				US0468-RFPKGHR
				US0468-RPKG
				US0468-RPKGHR
				US0468PL-RFPKG



Image	SKU	Price	Description	Item No. 2
-------	-----	-------	-------------	------------



US0468PL-RFPKGHR



US0468PL-RPKG



US0468PL-RPKGHR



US0458-220V

PRODUCT DESCRIPTION

Patriot Homecare Beds Full-Electric/Low Beds



Introducing the new Patriot Full-Electric Homecare Bed with a grid sleep surface!

Designed and engineered by the same team who develop the GF Long Term Care beds.

Sleep surface is MADE IN THE USA!

The single-motor design homecare bed is now available in a full-electric/low version. The Patriot Full-Electric Homecare Bed offers the same features as the Patriot Semi-Electric Homecare Bed with the additional feature of electric height adjustment. This makes it easier for caregivers to raise and lower the bed for patient care and linen changes. The electric height adjustment also enables the patient to easily position the bed for maximum comfort. The Full electric model can be adjusted to be a conventional height homecare bed OR a low bed for better inventory management!

Total assembled bed weighs only 178 pounds (bed and motor).

Powerful DC motor for effortless height adjustment.

Hi-Lo motor comes pre-installed on foot section for easy set-up and disassembly.

All functions fully-integrated into convenient handheld pendant.

Welded motor-guard protects motor and universal joint from damage during shipment and deliveries.

Available with Fiberboard Decorative Walnut Bed-ends or High-impact Plastic Bed-ends.

Designed to meet the requirements of applicable HCPCS codes.

Maximum patient weight 350 lbs.

450 lb evenly distributed total weight capacity including patient, mattress, rails, and accessories.

Height Adjustment:

Casters in Low Position:

Min. 9.5", Max. 20" (at deck level)

Casters in Standard Position:

Min. 13", Max. 23.5" (at deck level)

Overall Bed Size: 87" L x 36" W.

Lift Range: Head - 60°, Knee - 30°.

No external pinch points!

Slots in frame to ensure proper rail placement.

Locking pins secure bed ends to sleep deck. Attached lanyards prevent loss while moving or storing.

Only homecare bed to offer a 4-inch extension kit that can extend either the head, foot, or both head and foot decking. Bed length will be 84" with one kit, 88" with 2 kits.



The only homecare bed to offer optional corner mattress retainers.

Grid Deck offers the following advantages:

Much stronger frame.

Better mattress support.

No spring sag or bottoming out with heavier patients.

Easier to clean.

No connecting links or springs to lose.

NOTE:

Half-Rail A (GF6650A-1) will not work on the Full-Electric Patriot homecare bed when the bed is used as a Low Bed. You would be unable to lower the rail completely.

Half-Rail B (GF6580A-1) can be used on the Full-Electric Patriot homecare bed when used at conventional height or as a low-bed.

Warranty:

Two year warranty for defects in workmanship and materials of mechanical components, frame, and electronics.

During the warranty period, defective items will be repaired or replaced at manufacturer's option. The warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to the manufacturer.

Bed shipped in three cartons:

1. Foot section and motor.
2. Head section and Hi/Lo rod.
3. Bed end set.

US0458 Patriot Full-Electric Bed Only

US0458-RFPKG Patriot Full-Electric Bed with 1633-Reversible Foam Mattress and Full Chrome Rails

US0458-RPKG Patriot Full-Electric Bed with No Mattress and Full Chrome Rails

US0458PL Patriot Full-Electric Bed Only, Plastic Ends

US0458PL-220V Patriot Full-Electric Bed Only - 220 V - International Only, Plastic Ends

US0458PL-RFPKG Patriot Full-Electric Bed with 1633-Reversible Foam Mattress and Full Chrome Rails, Plastic Ends

US0458PL-RFPKGHA Patriot Full-Electric Bed with 1633-Reversible Foam Mattress and Clamp-On Half Chrome Rails, Plastic Ends

US0458PL-RPKG Patriot Full-Electric Bed with No Mattress and Full Chrome Rails, Plastic Ends

US0458PL-RPKGHA Patriot Full-Electric Bed with No Mattress and Clamp-On Half Chrome Rails, Plastic Ends

US0468PL-RFPKG Patriot Full-Electric Bed with 1633-Rolled Foam Mattress and and FDA



Compliant Full Rails, Plastic Ends

US0468PL-RFPKGHR Patriot Full-Electric Bed with 1633-Rolled Foam Mattress and and FDA Compliant Half No-Gap Rails, Plastic Ends

US0468PL-RPKG Patriot Full-Electric Bed with No Mattress and FDA Compliant Full Rails

US0468PL-RPKGHR Patriot Full-Electric Bed with No Mattress and FDA Compliant Half No-Gap Rails

US0458-220V Patriot Full-Electric Bed Only - 220 V - International Only

US0468-PKG Full-Electric Bed with 1633-Innerspring Mattress and FDA Full Rails

US0468-RFPKG Full-Electric Bed with 1633- Rolled Foam Mattress and FDA Full Rails

US0468-RFPKGHR Full-Electric Bed with 1633- Rolled Foam Mattress and FDA Half Rails

US0468-RPKG Full-Electric Bed with No Mattress and FDA Full Rails

US0468-PL-RPKGHR Full-Electric Bed with No Mattress and FDA Half Rails

ADDITIONAL INFORMATION

Item No. 2

1719, 1737, 1720, 1721, 1722, 1723, 1724, 1725, 1726, 1727,
1728, 1729, 1730, 1731, 1732, 1733, 1734, 1735, 1736