



## PATRIOT LX FULL-ELECTRIC HOMECARE BED

Patriot LX Full-Electric Homecare Bed

Technical Specifications

Total bed weight assembled: 126 lb

Maximum patient weight 350 lbs EVENLY DISTRIBUTED

450 lb distributed total weight capacity

including patient, mattress, rails and accessories

Overall Bed Size: 84" L x 35" W (80" x 35" sleep surface)

Lift Range: Head - 60, Knee - 30

Bed Height: Min. 10", Max. 20" (at deck level)



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [BED](#), [Full Electric Beds](#)

**Tags:** [full-electric bed](#), [Homecare bed](#), [patient bed](#)

### ADDITIONAL INFORMATION

<b>Item No. 2</b>	1709, 1701, 1702, 1703, 1765, 1704, 1705, 1706, 1707, 1708, 1766
-------------------	---

### VARIATIONS



Image	SKU	Price	Description	Item No. 2
				US6000-SFPKGQR
				US6000-RFPKGQR
				US6000
				US6000-RFPKG



Image	SKU	Price	Description	Item No. 2
				US6000-RFPKGHR
				US6000-RPKG
				US6000-RPKGHR
				US6000-RPKGQR



Image	SKU	Price	Description	Item No. 2
				US6000-SFPKG
				US6000-SFPKGHR
				690-8084-000

## PRODUCT DESCRIPTION

### Patriot LX Full-Electric Homecare Bed

Robust steel grid deck sleep surface provides superior mattress support and easy to clean surface compared to traditional spring decks

Head and foot sections can be raised or lowered simultaneously or independently with one hand control

Recessed main frame offers maximum pinch point protection between grid deck and frame - reducing the risk of injury

Lightweight design weighs only 126 lbs - lightest bed on the market

Bed setup is so simple that it can save on setup time and requires minimal maintenance



keeping your technicians time efficient

Patent pending integrated HI/LO mechanism makes bed height adjustment easy and reduces potential loss of parts such as HI/LO rod when used as rentals

Features dual HI/LO motors to adjust bed height from a low 10" to 20" easy and even allows for trendelenburg and slant positioning

Patent pending rail slots in frame ensure quick and proper rail placement

Single motor unit is designed for easy installation, requires less maintenance and provides the ultimate in quiet, smooth head and foot articulation

Attractive, "home-like" thermo fused bed ends offer a less institutional look and eliminate a bed component that often requires frequent repair and maintenance

Features SilverSolutions- an anti-microbial technology built into the paint that protects against bacteria and odor

During a power failure, the motor can lower both the head and foot section with just a single 9 volt battery

4" extension kit can be installed on head, foot or both to extend bed length to 84" or 88"

Proudly made in the U.S.A.

Grid Deck offers the following advantages:

Much stronger frame

Better mattress support

No spring sag or bottoming out with heavier patients

Easier to clean

No connecting links or springs to lose

Warranty:

Bed Frame: 5 Years

Electrical Components: 2 Years

All Other Durable Components: 1 Year

Technical Specifications

Total bed weight assembled: 126 lb

Maximum patient weight 350 lbs EVENLY DISTRIBUTED

450 lb distributed total weight capacity

including patient, mattress, rails and accessories

Overall Bed Size: 84" L x 35" W (80" x 35" sleep surface)

Lift Range: Head - 60, Knee - 30

Bed Height: Min. 10", Max. 20" (at deck level)

US6000 Patriot LX Full-Electric Bed Only

US6000-RFPKG Patriot LX Full-Electric Bed with 1633-Rolled Foam Mattress and FDA Full



Rails  
US6000-RFPKGHR Patriot LX Full-Electric Bed with 1633-Rolled Foam Mattress and FDA Half  
Rails  
US6000-RFPKGQR Patriot LX Full-Electric Bed with 1633-Rolled Foam Mattress and FDA  
Quarter Rails  
US6000-RPKG Patriot LX Full-Electric Bed with No Mattress and FDA Full Rails  
US6000-RPKGHR Patriot LX Full-Electric Bed with No Mattress and FDA Half Rails  
US6000-RPKGQR Patriot LX Full-Electric Bed with No Mattress and FDA Quarter Rails  
US6000-SFPKG Patriot LX Full-Electric Bed with 1633-Standard Foam Mattress & FDA Full  
Rails  
US6000-SFPKGHR Patriot LX Full-Electric Bed with 1633-Standard Foam Mattress & FDA Half  
Rails  
US6000-SFPKGQR Patriot LX Full-Electric Bed with 1633-Standard Foam Mattress & FDA  
Quarter Rails  
690-8084-000 Extension Kit for Patriot LX

## **ADDITIONAL INFORMATION**

<b>Item No. 2</b>	1709, 1701, 1702, 1703, 1765, 1704, 1705, 1706, 1707, 1708, 1766
-------------------	---