



## PEDIATRIC STETHOSCOPE

Pediatric Stethoscope

507 Black

507LB Lt. Blue

507P Pink



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [DIAGNOSTIC INSTRUMENTS](#), [Stethoscopes](#)

**Tags:** [pediatric stethoscope](#), [stethoscope](#)



### ADDITIONAL INFORMATION

<b>Weight</b>	N/A
<b>Dimensions</b>	N/A
<b>Color</b>	107, 257, 242

### VARIATIONS



**Image**

**SKU**

**Price**

**Description**

**Color**



Black



Slate Blue



Pink

**GALLERY IMAGES**





## PRODUCT DESCRIPTION

Pediatric Stethoscope

507 STETH PED DUAL HEAD BLK Black

507LB STETH PED DUAL HEAD LT BLUE Lt. Blue

507P STETH PED DUAL HEAD PINK Pink

Lightweight anodized aluminum rotating dual head chestpiece with recessed non-chill ring and adjustable chrome-plated brass external spring binaural. 22" latex-free tubing.

## ADDITIONAL INFORMATION

<b>Weight</b>	N/A
<b>Dimensions</b>	N/A
<b>Color</b>	107, 257, 242