



5711 PILL CUTTER

5711 Pill Cutter

Splits tablets up to 9/16" diameter without crumbling.
Features a storage compartment to hold unused portions.



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Daily Living Aids](#), [Personal Care/Daily Living Aids](#)

Tags: [pill cutter](#), [pill splitter](#)

1

PRODUCT DESCRIPTION

5711 Pill Cutter

Sturdy plastic cutter with stainless steel edge requires little pressure to split coated or uncoated pills, caplets and tablets.

Splits tablets up to 9/16" diameter without crumbling.

Features a storage compartment to hold unused portions.

1