




## PLATINUM COLLECTION BATH SEATS

Platinum Collection Bath Seats  
7921A-4 Standard Grey  
7923A Platinum Grey  
7931A Standard Grey



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [Bath Seats](#), [BATHROOM SAFETY](#)


**Tags:** [collection bath seat](#), [platinum bath seat](#)

### ADDITIONAL INFORMATION

<b>Weight</b>	N/A
<b>Dimensions</b>	N/A
<b>Item No. 2</b>	558, 559, 560
<b>Color</b>	<a href="#">Grey</a> , <a href="#">Platinum Grey</a>

### VARIATIONS



Image	SKU	Price	Description	Item No. 2
				7921A-4
				7923A
				7931A

**PRODUCT DESCRIPTION**

Platinum Collection Bath Seats

7921A-4 Lumex Platinum Collection Bath Seat with Backrest Standard Grey 4 Per Case

7923A Lumex Platinum Collection Bath Seat With Backrest and Seat Cutout Platinum Grey

7931A Lumex Platinum Collection Bath Seat without Backrest Standard Grey

Anodized aluminum frame is lightweight, durable and rust-resistant

Seat height is adjustable in 1" increments

Attractive platinum grey color

Available with or without backrest

Comes assembled with tool-free back attachment



Includes laminated color operating instructions  
Designed to meet the requirements of HCPCS code: E0245  
Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED  
Lifetime limited warranty  
Specifications  
Overall Width 20"  
Overall Depth 17"  
Seat Width 20"  
Seat Depth 12"  
Minimum Seat Height 14"  
Maximum Seat Height 21"  
Maximum Weight Capacity 350 lb EVENLY DISTRIBUTED

#### **ADDITIONAL INFORMATION**

<b>Weight</b>	N/A
<b>Dimensions</b>	N/A
<b>Item No. 2</b>	558, 559, 560
<b>Color</b>	<a href="#">Grey</a> , <a href="#">Platinum Grey</a>