




POLYPROPYLENE ISOLATION GOWN 10 / BAG

Polypropylene Isolation Gown
10 / BAG
Water / Liquid resistant
Latex Free



SKU: N/A

Price:

Stock: 15 in stock (can be backordered)

Categories: [Isolation Gowns](#), [PPE](#)

Tags: [apparel](#), [gown](#)

1



PRODUCT DESCRIPTION



Polypropylene Isolation Gown

10 / BAG

- Isolation gown
- Non-woven polypropylene in yellow
- Water / Liquid resistant
- 30 gm/m²
- 127 cm x 152 cm
- Latex Free