



## HOM549 PORTABLE DIGITAL BABY/ SMALL ANIMAL SCALE

HOM549 Portable Digital Baby/ Small Animal Scale  
44lb/20 kg weight capacity.  
Operated by 1 9V battery (included)



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#), [Scales](#)

**Tags:** [digital animal scale](#), [digital baby scale](#)

1

## PRODUCT DESCRIPTION

HOM549 Portable Digital Baby/ Small Animal Scale

Portable digital pediatric scale with motion sensing weighing technology for faster readings of active infants. Features a 6/10" LCD display showing weight in pounds and fractions of ounces or in kilograms.

44lb/20 kg weight capacity.

Includes tare function to remove pad or blanket weight.

Operated by 1 9V battery (included)

1