




RECLINING SHOWER CHAIR

Reclining Shower Chair

89330 Lumex 20" PVC Reclining Shower Commode Chair Specify Color

89440 Lumex 26" Bariatric Reclining Commode Bath Seat Specify Color

89550 Lumex 30" Bariatric Reclining Commode Bath Seat Specify Color



SKU: N/A

Price:

Stock: Available on backorder

Categories: [BATHROOM SAFETY](#), [Shower Transport Chairs](#)




Tags: [reclining chair](#), [shower chair](#)

ADDITIONAL INFORMATION

Weight	6 lbs
Dimensions	50 × 40 × 50 in
Item No. 4	562, 563, 564

VARIATIONS



Image	SKU	Price	Description	Item No. 4
				89330
				89440
				89550

PRODUCT DESCRIPTION

Reclining Shower Chair

89330 Lumex 20" PVC Reclining Shower Commode Chair with Footrest, Swing Away Arms, Soft Seat and Safety Belt Specify Color

89440 Lumex 26" Bariatric Reclining Commode Bath Seat with Elevated Leg Rest, Swing Away Arms, Soft Seat, and Safety Belt Specify Color

89550 Lumex 30" Bariatric Reclining Commode Bath Seat with Elevated Leg Rest, Swing Away Arms, Soft Seat, and Safety Belt Specify Color



Constructed of the highest quality health care grade PVC
Threaded stem casters for additional safety.
Ergonomic curved rails provide support without sharp edges and reclining chairs have anti-slip handgrips for added safety.
One year limited warranty on PVC frame; six month limited warranty on all other materials.
Deluxe elongated extra heavy-duty, contoured open front soft seat.
Padded lumbar support and head section.
Dual swing away arms for patient transfer can be secured with red safety straps
Elevated Leg Rest (89440 and 89550 only)
Mesh available in your choice of 3 colors -
Royal Blue mesh will be provided for online orders. If you desire another color, please order via Customer Service at 1-800-347-56781-800-347-5678.

ADDITIONAL INFORMATION

Weight	6 lbs
Dimensions	50 × 40 × 50 in
Item No. 4	562, 563, 564