



SAFESHIELD™ DISPOSABLE STERILE SCALPELS 10/BOX

Safeshield Disposable Sterile Scalpels

10/BOX

2980#10 Size 10 Blue Handle with Brown Tab

2980#11 Size 11 Blue Handle with Pink Tab

2980#14 Size 14 Blue Handle with Red Tab

2980#15 Size 15 Blue Handle with Purple Tab


2980#15C Size 15C Blue Handle with Blue Tab

2980#20 Size 20 Blue Handle with Light Brown Tab

2980#21 Size 21 Blue Handle with Green Tab

2980#23 Size 23 Blue Handle with Turquoise Tab

2980#24 Size 24 Blue Handle with Grey Tab



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Instruments](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tags: [disposable scalpels](#), [sterile scalpels](#)

ADDITIONAL INFORMATION

Weight 4.50 lbs

Dimensions 4.80 × 3.00 × 2.00 in

Size #11, #14



PRODUCT DESCRIPTION

Safeshield Disposable Sterile Scalpels

2980#10 Feather® Safeshield™ Disposable Sterile Scalpels, Size 10 Blue Handle with Brown Tab

2980#11 Feather® Safeshield™ Disposable Sterile Scalpels, Size 11 Blue Handle with Pink Tab

2980#14 Feather® Safeshield™ Disposable Sterile Scalpels, Size 14 Blue Handle with Red Tab

2980#15 Feather® Safeshield™ Disposable Sterile Scalpels, Size 15 Blue Handle with Purple Tab

2980#15C Feather® Safeshield™ Disposable Sterile Scalpels, Size 15C--Specialty Blade (Plastics) Blue Handle with Blue Tab

2980#20 Feather® Safeshield™ Disposable Sterile Scalpels, Size 20 Blue Handle with Light Brown Tab

2980#21 Feather® Safeshield™ Disposable Sterile Scalpels, Size 21 Blue Handle with Green Tab

2980#23 Feather® Safeshield™ Disposable Sterile Scalpels, Size 23 Blue Handle with Turquoise Tab

2980#24 Feather® Safeshield™ Disposable Sterile Scalpels, Size 24 Blue Handle with Grey Tab

Feather® Safeshield™ Disposable Sterile Scalpels are ideal for operating rooms and surgical facilities and are excellent for delicate procedures when precise control is needed. Feather, a blade manufacturer of global renown, utilizes high-precision technology and carefully selected materials. Each Feather scalpel is made in Japan and meets strict ISO standards, resulting in exceptionally high quality and consistent strength and sharpness. Designed with an ergonomically rounded SUREgrip™ plastic handle that fits the hand comfortably and outfitted with a stainless steel blade, Safeshield™ Safety Scalpels feature a locking retractable protective shield that helps prevent injuries before, during, and after use.

Features:

- Sterile: Gamma-sterilized and ready for use immediately after opening
- Single Use/Disposable: Helps prevent patient cross contamination
- Individually packaged: Equipped with clear window for product selection verification and pull tabs for easy opening



- Clear locking retractable safety shield: Can be operated with one hand
- Available in 11 sizes: Including most common blade sizes and select specialty blade size
- Color-coded by size: Release tab on protective shield and coordinating packaging assists with blade identification
- Ruler printed on handle: Includes cm and mm markings
- Similar weight as conventional disposable scalpel
- Not made with natural rubber latex
- ISO 9001 and 13485 certified; ISO 14001 compliant

Note: Dispose of scalpel in an approved sharps container in accordance with your facility's sharps disposal procedures

ADDITIONAL INFORMATION

Weight	4.50 lbs
Dimensions	4.80 × 3.00 × 2.00 in
Size	#11, #14