



## SPONGE BOWLS


Sponge Bowls

3248 Capacity: 12 oz.

3249 Capacity: 1 qt.

3250 Capacity: 1 3/8 qts.

3251 Capacity: 1 7/8 qts.



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#),  
[Stainless Steelware](#)

**Tags:** [medical bowls](#), [sponge bowls](#)

## ADDITIONAL INFORMATION

<b>Weight</b>	N/A
<b>Dimensions</b>	N/A
<b>Item No. 4</b>	234, 235, 236, 237

## VARIATIONS



Image	SKU	Price	Description	Item No. 4
				3248
				3249
				3250
				3251

**PRODUCT DESCRIPTION**

Sponge Bowls

3248 4 1/2" (width) x 2" (Depth), Capacity: 12 oz.



3249 5 3/4" x 2 1/2". Capacity: 1 qt.  
3250 6 1/2" x 3". Capacity: 1 3/8 qts.  
3251 7 3/8" x 3". Capacity: 1 7/8 qts.  
18-8 (304) Stainless Steel with seamless construction.

#### **ADDITIONAL INFORMATION**

<b>Weight</b>	N/A
<b>Dimensions</b>	N/A
<b>Item No. 4</b>	234, 235, 236, 237