



THERMOMETER, DISPOSABLE, 60 SECOND, BULK

Thermometer, Disposable, 60 Second, Bulk
1858 60 SECOND, BULK 24EA/BX Rigid Tip
1858CASE 60 SECOND, BULK, 240/CS RIGID TIP



SKU: N/A

Price:

Stock: Out of stock

Categories: [MEDICAL, SURGICAL and PHYSICIAN, Thermometers](#)



Tags: [disposable thermometer, thermometer](#)

ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Item No. 2	566, 565

VARIATIONS



Image	SKU	Price	Description	Item No. 2
				1858CASE
				1858

PRODUCT DESCRIPTION

Thermometer, Disposable, 60 Second, Bulk

1858 DIGITAL THERMOMETER, 60 SECOND, BULK 24EA/BX Rigid Tip

1858CASE DIGITAL THERMOMETER, 60 SECOND, BULK, 240/CS RIGID TIP

Specifications:

Range: 90.0°F-109.9°F

Accuracy: +/- 0.2°F, 95.9°F-107.6°F

+/- 0.4°F under 95.9°F or over 107.6°F

Memory: Last measured Value

Battery Life: Approximately 200 hours

Dimension: 12.4cm (L) x 1.8cm (W) x 0.9cm (H)

Weight: Approximately 10 grams including battery

ADDITIONAL INFORMATION



Weight	N/A
Dimensions	N/A
Item No. 2	566, 565