



1325 TUNING FORK – ADJUSTABLE FREQUENCY

1325 Tuning Fork - Adjustable Frequency
Tines 3/8", Length 8"



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Instruments](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tags: [adjustabl tuning fork](#), [tuning fork](#)

1

PRODUCT DESCRIPTION

1325 Tuning Fork - Adjustable Frequency
Tines 3/8", Length 8"

All aluminum tuning fork, with adjustable weight scale in seven frequencies: C128, D144, E160, F170.7, G192, A213.3 and B240 (weights removed). Scale letters are stamped on stems.

1