



1325-1 TUNING FORK – SENSORY EVALUATION

1325-1 Tuning Fork – Sensory Evaluation

Tines 3/8", Handle 2 7/8", Length 11 1/2"



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Instruments](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tags: [evaluation fork](#), [tuning fork](#)

1

PRODUCT DESCRIPTION

1325-1 Tuning Fork – Sensory Evaluation

Tines 3/8", Handle 2 7/8", Length 11 1/2"

For use on damaged regenerating nerves. Aluminum alloy 30 Hz tuning fork with massive steel weights. Machined from 1 1/4" bar stock.

1