



TUNING FORKS

Tuning Forks

1313 Set of all 5


1314 C128 fixed weight

1315 C256 fixed weight

1316 C512 without weight

1317 C1024 without weight

1318 C2048 without weight



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Instruments](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tag: [tuning forks](#)

ADDITIONAL INFORMATION

Item No. 1 148, 149, 150, 151, 152, 153

VARIATIONS







Image	SKU	Price	Description	Item No. 1
				1313
				1314
				1315
				1316



Image	SKU	Price	Description	Item No. 1
				1317
				1318

PRODUCT DESCRIPTION

Tuning Forks

1313 Set of all 5 Comes with cushioned vinyl zippered case

1314 C128 fixed weight

1315 C256 fixed weight

1316 C512 without weight

1317 C1024 without weight

1318 C2048 without weight

Aluminum alloy tuning forks: corrosion and tarnishing resistant

Tines are 1¾" long and 3/8" diameter

Single frequency pure tone will last long

ADDITIONAL INFORMATION

Item No. 1 148, 149, 150, 151, 152, 153

