



VACUTAINER NEEDLE HOLDER SINGLE-USE

Vacutainer Needle Holder Single-Use

364815 BD Vacutainer one-use holder 250 (per bag)



SKU: N/A

Price:

Stock: Out of stock

Category: [PPE](#)

1



PRODUCT DESCRIPTION

Vacutainer Needle Holder Single-Use

Clear Vacutainer Needle Holders for single-use are standard plastic holders that accept disposable and reusable double-ended needles with threaded hubs. Holders convert BD Vacutainer tube and needle into a syringe-like assembly that facilitates use. Not autoclavable.



Helps customers to comply with the OSHA safety and health information bulletin on the disposal of contaminated needles and blood tube holders. Non-Sterile.

364815 BD Vacutainer one-use holder 250 (per bag)