




VACUTAINER NEEDLE HOLDER SINGLE USE

Vacutainer Needle Holder Single Use

Clear Vacutainer Needle Holders for single use are standard plastic holders that accept disposable and reusable double-ended needles with threaded hubs.



SKU: N/A

Price:

Stock: 1000 in stock (can be backordered)

Categories: [Accessories](#), [IV THERAPY](#)



ADDITIONAL INFORMATION



Weight N/A

Dimensions N/A

Package Pack of 50, Pack of 100, Bag of 250

VARIATIONS



Image	SKU Price	Description	Package
		<p>Vacutainer Needle Holder Single Use Pack of 100 Clear Vacutainer Needle Holders for single use are standard plastic holders that accept disposable and reusable double-ended needles with threaded hubs.</p>	Pack of 100
		<p>Vacutainer Needle Holder Single Use Pack of 250 Clear Vacutainer Needle Holders for single use are standard plastic holders that accept disposable and reusable double-ended needles with threaded hubs.</p>	Bag of 250

PRODUCT DESCRIPTION

Vacutainer Needle Holder Single Use

Clear Vacutainer Needle Holders for single use are standard plastic holders that accept disposable and reusable double-ended needles with threaded hubs. Holders convert BD Vacutainer tube and needle into a syringe-like assembly that facilitates use. Not autoclavable. Helps customers to comply with the OSHA safety and health information bulletin on the disposal of contaminated needles and blood tube holders. Non-Sterile.

ADDITIONAL INFORMATION

Weight	N/A
Dimensions	N/A
Package	Pack of 50, Pack of 100, Bag of 250